

Set of four mattress side retainers – fitting guide

		
<p>Included are 4 side retainers, 16 screws and fitting guide.</p>	<p>Each side retainer is two components.</p>	<p>These have already been clipped together for your convenience.</p>
		
<p>Position two retainers each side, one near the shoulder area and one by the thigh</p>	<p>Position the retainer on the inner frame at a point that will hold the mattress as it is adjusted</p>	<p>Fasten the two screws into the frame and push the face plate flush to the face of the bed base.</p>
		
<p>Reposition the mattress. Pull out or push in the faceplates to hold the mattress nice and snug.</p>	<p>Then fasten the two parts of the retainer securely using the remaining eight screws.</p>	